


Single-use
without a Single
Compromise



GlideScope Spectrum QC

New generation Spectrum™ QC video laryngoscopes couple innovative magnetic QuickConnect™ technology with imaging improvements that empower clinicians to act with confidence in routine and emergent situations.



Innovative Connections

Magnetic QuickConnect for plug-and-play simplicity

Fast, Efficient Care

QuickConnect for Spectrum QC and BFlex allows you to unify your GlideScope® Core™ system with one innovative video cable solution for procedural efficiency



New connector, same broad portfolio of nine blade sizes and types.
For more information visit verathon.com/spectrumqc



Enhanced Airway Views

50% improvement in airway image sharpness and vividness¹











Expanded Visualization

Latest design expands visualization in the airway with 14% greater field of view²

1. Comparing Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum HDMI LoPro/Mac 3 & 4
2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac

Spectrum Product Comparison

Workflow and Connectivity	Spectrum™ QC (New)	Spectrum™
Connector Type	Magnetic QuickConnect™ 	HDMI-like 
Connector Retention Force	50% higher retention force ¹	–
Dual-View Cable Requirements	Two GlideScope® Core™ QuickConnect Cables	One Core Smart Cable One Core QuickConnect Cable
Imaging		
Field of View	14% greater field of view ²	–
Dynamic Light Control®	Yes	Yes
LED Brightness	Mac 30% brighter; Miller 100% brighter ³	–
Image Sharpness and Vividness	50% sharper / more-vivid ⁴	–
Visualization System Compatibility		
GlideScope Core Monitors 	GlideScope Core QuickConnect Cable (for Spectrum QC and BFlex) 	GlideScope Core Smart Cable (for Spectrum) 
GlideScope Video Monitor (GVM) 	Not compatible	Spectrum Smart Cable 
GlideScope® Go™ 	GlideScope Go 2 (New)	GlideScope Go

1. Spectrum QC with Core QuickConnect cable versus Spectrum with Core Smart Cable
 2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac
 3. Comparing Spectrum QC Mac/Miller versus Spectrum Mac/Miller
 4. Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum LoPro/Mac 3 & 4

To learn more, contact your local representative.



verathon.com

Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA
Tel: +1 800 331 2313
(USA & Canada only)
Tel: +1 425 867 1348
Fax: +1 425 883 2896



Verathon Medical (Europe) B.V.
Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia
Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970
Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201